

The Waterpik® Water Flosser: Significantly more effective than interdental brush for improving gingival health!

Water flosser compared to interdental brush on bleeding scores and gingival abrasion

Slot DE, Lyle DM, Van der Sluijs E, Hennequin-Hoenderdos N, Van der Weijden F. J Dent Res 2018; 97(Special Iss. B): Abstract #0622 (www.iadr.org). Conducted at Academic Center for Dentistry Amsterdam (ACTA), Netherlands.

Objective

To compare the effectiveness of a water flosser (WF) and interdental brush (IDB) on bleeding indices and gingival abrasion.

Methodology

Seventy-eight subjects completed this 4-week, randomised controlled trial. Subjects were assigned to one of two groups; Waterpik® Water Flosser (WF) plus a manual toothbrush or interdental brush (IDB) plus a manual toothbrush. Gingival inflammation was evaluated by measuring Bleeding on Pocket Probing (BOPP) and Bleeding on Marginal Probing (BOMP). Data was collected on contra-lateral quadrants. The Gingival Abrasion Score (GAS) was used to compare the incidence of abrasion between the groups.


Results

Both groups demonstrated a significant reduction in BOPP and BOMP from baseline to 4 weeks for all sites and interdental sites separately. The WF group was significantly more effective than the IDB group for reducing BOPP for all sites at W4 ($p=0.030$) and BOMP for all sites and interdental sites at W4 ($p=0.003$, $p=0.019$ respectively). There were no differences in gingival abrasion scores between the groups.

Conclusion

The Waterpik® Water Flosser is significantly more effective than the interdental brush for improving gingival health in this clinical study.

Bleeding on Pocket Probing
(4-week results for all sites)



Bleeding on Marginal Probing
(4-week results for all sites)

