

mis[®]
CONICAL CONNECTION
PROSTHETIC OPTIONS

MIS conical connection implants are the outcome of an exceptionally high-level R&D process resulting in implant systems that are simple, easy-to-use with many enhanced functionality and performance features. the conical connection concave abutment profile helps gain thicker soft tissue for better esthetic results. In addition, concave abutments have a special groove to verify proper positioning.







V3
V3-10330
V3-11330
V3-13330
V3-16330

Cover Screw & Healing Caps



CC1-00277



CN-HS340 (Ø4mm H=3mm)
CN-HS440 (Ø4mm H=4mm)
CN-HS540 (Ø4mm H=5mm)
CN-HS640 (Ø4mm H=6mm)
CN-HS840 (Ø4mm H=8mm)



CN-HA248 (Ø4.8mm H=2mm)
CN-HA348 (Ø4.8mm H=3mm)
CN-HA448 (Ø4.8mm H=4mm)
CN-HA548 (Ø4.8mm H=5mm)
CN-HA648 (Ø4.8mm H=6mm)
CN-HA848 (Ø4.8mm H=8mm)

Impression Copings



VN-IO331 (Ø3.3mm H=11mm)
VN-IO335 (Ø3.3mm H=15mm)



VN-IO401 (Ø4mm H=11mm)
VN-IO405 (Ø4mm H=15mm)



VN-IC401 (Ø4mm H=11mm)
VN-IC405 (Ø4mm H=15mm)



VN-IC400 (Ø4mm)




VN-PF340 (Ø4mm)

Analog



VN-RSM10


Temporary Cylinders



VN-TPH50 (Ø3.5mm)



VN-TPA10 (Ø4mm H=1.5mm)



CN-TC013 (Ø4mm H=1mm)
CN-TC023 (Ø4mm H=2mm)
CN-TC033 (Ø4mm H=3mm)





VN-TCl13 (Ø4mm H=1mm)
VN-TCl23 (Ø4mm H=2mm)
VN-TCl33 (Ø4mm H=3mm)




Fixed Prosthetic Options

Prosthetic Options

Gold Plastic Cylinders





Abutments for Cement Retained Restorations





* Includes burnout anti-rotation plastic cap MN-ICH40

Angulated Abutment



Zirconia Solution



*All superstructures are provided with the relevant prosthetic screws

NARROW

Multi-Unit



CM-N1480 (H=1mm)
CM-N2480 (H=2mm)
CM-N3480 (H=3mm)
CM-N4480 (H=4mm)
CM-N5480 (H=5mm)

Ball Attachments



CN-B1330 (H=1mm)
CN-B2330 (H=2mm)
CN-B3330 (H=3mm)
CN-B4330 (H=4mm)

OT-Equators



CK-NOE1 (H=1mm)
CK-NOE2 (H=2mm)
CK-NOE3 (H=3mm)
CK-NOE4 (H=4mm)
CK-NOE5 (H=5mm)

* Available in kits only

Ti-Base



VN-TB001 (Ø4.4mm H=0.5mm)
VN-TB015 (Ø4.4mm H=1.5mm)



CN-TB002 (Ø4.2mm H=0.5mm)
CN-TBF15 (Ø4.2mm H=1.5mm)



VN-IN060 (Ø4.4mm H=0.5mm)
VN-IN160 (Ø4.4mm H=1.5mm)




VN-RSM10



GI-SCH11

Implant Details

NARROW



3.30
Diameter

*All superstructures are provided
with the relevant prosthetic screws

- ☐ Free rotation
- ☐ Anti-rotation



C1
C1-10330
C1-11330
C1-13330
C1-16330

Fixed Prosthetic Options

Cover Screw & Healing Caps



CN1-00330



CN-HS340 (Ø4mm H=3mm)
CN-HS440 (Ø4mm H=4mm)
CN-HS540 (Ø4mm H=5mm)
CN-HS640 (Ø4mm H=6mm)
CN-HS840 (Ø4mm H=8mm)



CN-HA248 (Ø4.8mm H=2mm)
CN-HA348 (Ø4.8mm H=3mm)
CN-HA448 (Ø4.8mm H=4mm)
CN-HA548 (Ø4.8mm H=5mm)
CN-HA648 (Ø4.8mm H=6mm)
CN-HA848 (Ø4.8mm H=8mm)

Impression Copings



CN-IO331 (Ø3.3mm H=11mm)
CN-IO335 (Ø3.3mm H=15mm)



CN-IO330 (Ø4mm H=11mm)
CN-IO405 (Ø4mm H=15mm)



CN-IC401 (Ø4mm H=11mm)
CN-IC405 (Ø4mm H=15mm)



CN-IT300 (Ø4mm)




CN-PF330

Analog



CN-RSM10

Temporary Cylinders



CN-TPH50 (Ø3.5mm)



CN-TPA10 (Ø4mm H=1.5mm)



CN-TC013 (Ø4mm H=1mm)
CN-TC023 (Ø4mm H=2mm)
CN-TC033 (Ø4mm H=3mm)



CN-TC113 (Ø4mm H=1mm)
CN-TC123 (Ø4mm H=2mm)
CN-TC133 (Ø4mm H=3mm)

Prosthetic Options

Prosthetic screw: Standard



MN-S0160

Gold Plastic Cylinders



CN-GPC10 (Ø4.4mm H=1mm)



CN-GP010 (Ø4.4mm H=1mm)

Abutments for Cement Retained Restorations



CN-MAC10 (Ø4mm H=1mm)
CN-MAC30 (Ø4mm H=3mm)



CN-AN001 (Ø4.8mm H=1mm)
CN-AT002 (Ø4.8mm H=2mm)
CN-AN003 (Ø4.8mm H=3mm)



CN-CPK41 (C=4mm H=1mm)
CN-CPK42 (C=4mm H=2mm)
CN-CPK43 (C=4mm H=3mm)
CN-CPK44 (C=4mm H=4mm)

CN-CPK61 (C=6mm H=1mm)
CN-CPK62 (C=6mm H=2mm)
CN-CPK63 (C=6mm H=3mm)
CN-CPK64 (C=6mm H=4mm)



MN-IC040



MN-ICH40

* Includes burnout anti-rotation plastic cap MN-ICH40

Angulated Abutment




CN-AN010 (Ø4.8mm H=1mm)
CN-AT102 (Ø4.8mm H=2mm)
CN-AN310 (Ø4.8mm H=3mm)



CN-AN020 (Ø4.8mm H=1mm)
CN-AT202 (Ø4.8mm H=2mm)
CN-AN320 (Ø4.8mm H=3mm)

Zirconia Solution



Non-Cemented

CK-NIZ6 (H=0.5mm)



Pre-Cemented

CK-NIC6 (H=0.5mm)




*All superstructures are provided with the relevant prosthetic screws

NARROW

Multi-Unit	Ball Attachments	OT-Equators	Ti-Base	
 CM-N1480 (Ø4.8mm H=1mm) CM-N2480 (Ø4.8mm H=2mm) CM-N3480 (Ø4.8mm H=3mm) CM-N4480 (Ø4.8mm H=4mm) CM-N5480 (Ø4.8mm H=5mm)	 CN-B1330 (H=1mm) CN-B2330 (H=2mm) CN-B3330 (H=3mm) CN-B4330 (H=4mm)	 CK-NOE1 (H=1mm) CK-NOE2 (H=2mm) CK-NOE3 (H=3mm) CK-NOE4 (H=4mm) CK-NOE5 (H=5mm) * Available in kits only	 CN-TB001 (Ø4.4mm H=0.5mm) CN-TB015 (Ø4.4mm H=1.5mm)	 CN-TB002 (Ø4.4mm H=0.5mm) CN-TBF15 (Ø4.4mm H=1.5mm)
			 CN-IN060 (Ø4.4mm H=0.5mm) CN-IN160 (Ø4.4mm H=1.5mm)	
			 CN-RSM10	 GI-SCH11

Implant Details

NARROW



3.30
Diameter



*All superstructures are provided
with the relevant prosthetic screws

- ☐ Free rotation
- ☒ Anti-rotation

STANDARD C1 (Ø3.75, Ø4.20) V3 (Ø3.90, Ø4.30, Ø5)

Restorative Options for Standard Platform

Cover Screw





Fixed Prosthetic Options



Healing Caps






Impression Copings



Analogs



Temporary Cylinders



Prosthetic Options

All superstructures are provided with the relevant prosthetic screws:

Standard



Optional gold



Gold Plastic Cylinders



CS-GPC10 (Ø4.4mm H=1mm)



CS-GP010 (Ø4.4mm H=1mm)

Abutments for Cement Retained Restorations



CS-A0010 (Ø4.8mm H=1mm)
CS-AT002 (Ø4.8mm H=2mm)
CS-A0030 (Ø4.8mm H=3mm)



VS-AT001 (Ø5.8mm H=1mm)
VS-AT002 (Ø5.8mm H=2mm)
VS-AT003 (Ø5.8mm H=3mm)



CS-MAC10 (Ø4.8mm H=1mm)
CS-MC103 (Ø4.8mm H=3mm)



VS-MC101 (Ø5.8mm H=1mm)
VS-MC103 (Ø5.8mm H=3mm)



CS-CPK41 (C=4mm H=1mm)
CS-CPK42 (C=4mm H=2mm)
CS-CPK43 (C=4mm H=3mm)
CS-CPK44 (C=4mm H=4mm)

CS-CPK61 (C=6mm H=1mm)
CS-CPK62 (C=6mm H=2mm)
CS-CPK63 (C=6mm H=3mm)
CS-CPK64 (C=6mm H=4mm)

CS-CPK81 (C=8mm H=1mm)
CS-CPK82 (C=8mm H=2mm)
CS-CPK83 (C=8mm H=3mm)
CS-CPK84 (C=8mm H=4mm)



VS-P6155 (C=6mm H=1mm)
VS-P6255 (C=6mm H=2mm)
VS-P6355 (C=6mm H=3mm)
VS-P6455 (C=6mm H=4mm)

VS-P8155 (C=8mm H=1mm)
VS-P8255 (C=8mm H=2mm)
VS-P8355 (C=8mm H=3mm)
VS-P8455 (C=8mm H=4mm)



MD-IC040



MD-ICH40



MW-IC040



MW-ICH40

* Includes burnout anti-rotation plastic cap MD-ICH40

Angulated Abutments



CS-A1510 (Ø4.8mm H=1mm)
CS-AT152 (Ø4.8mm H=2mm)
CS-A1530 (Ø4.8mm H=3mm)



CS-A2510 (Ø4.8mm H=1mm)
CS-AT252 (Ø4.8mm H=2mm)
CS-A2530 (Ø4.8mm H=3mm)



CS-AN151 (Ø5.8mm H=1mm)
VS-AT152 (Ø5.8mm H=2mm)
VS-AT153 (Ø5.8mm H=3mm)




CS-AN251 (Ø5.8mm H=1mm)
VS-AT252 (Ø5.8mm H=2mm)
VS-AT253 (Ø5.8mm H=3mm)

Zirconia Solution



Non-Cemented

CK-SIZ6 H=0.5mm



CK-SPZ4 H=0.5mm

* Available in kits only



Pre-Cemented

CK-SIC6 H=0.5mm



CK-SPC4 H=0.5mm

* Available in kits only

○ Free rotation

⬡ Anti-rotation

STANDARD

STANDARD


C1 (Ø3.75, Ø4.20)
V3 (Ø3.90, Ø4.30, Ø5)

Prosthetic Options (Multi-Unit, Ball Attachments, OT-Equators, Ti-Base)


Multi-Unit	Ball Attachments	OT-Equators	Ti-Base		
 CM-S1480 (H=1mm) CM-S2480 (H=2mm) CM-S3480 (H=3mm) CM-S4480 (H=4mm) CM-S5480 (H=5mm) Angulated  CM-SA171 CM-SA301	 CS-B1375 (H=1mm) CS-B2375 (H=2mm) CS-B3375 (H=3mm) CS-B4375 (H=4mm)	 CK-SOE1 (H=1mm) CK-SOE2 (H=2mm) CK-SOE3 (H=3mm) CK-SOE4 (H=4mm) CK-SOE5 (H=5mm) * Available in kits only	 CS-TB001 (H=0.5mm) CS-TB015 (H=1.5mm)	 CS-TB002 (H=0.5mm) CS-TBF15 (H=1.5mm)	
			 CS-IN060 (Ø4.4mm H=0.5mm) CS-IN160 (Ø4.4mm H=1.5mm)		
			Analog  CS-RSM10	Burnout  GI-SCH11	Scanning Abutment  CS-SP001

* Scanning Abutment is sold separately and is not part of the Ti-Base package.


Implant Details




3.75
Diameter




4.20
Diameter



3.90
Diameter



4.30
Diameter



5
Diameter

*All superstructures are provided
with the relevant prosthetic screws

- Free rotation
- ⬡ Anti-rotation



C1
C1-08500
C1-10500
C1-11500
C1-13500
C1-16500

Fixed Prosthetic Options

Cover Screw



CC1-00500

Healing Caps



CW-HS355 (Ø5.5, H=3mm)
CW-HS455 (Ø5.5, H=4mm)
CW-HS555 (Ø5.5, H=5mm)



CW-HA263 (Ø6.3, H=2mm)
CW-HA363 (Ø6.3, H=3mm)
CW-HA563 (Ø6.3, H=5mm)


Impression Copings



CW-IO501 (Ø5, H=11mm)
CW-IO505 (Ø5, H=15mm)



CW-IO500 (Ø5.5, H=11mm)
CW-IO555 (Ø5.5, H=15mm)



CW-IC551 (Ø5.5, H=11mm)
CW-IC555 (Ø5.5, H=15mm)



CW-IT300 (Ø5.5)



CW-PF500 (Ø5.5)

Analog



CW-RSM10

Prosthetic Options

Prosthetic screws:

Standard



Optional gold



Temporary Cylinder



CW-TC013 (Ø5.5, H=1mm)
CW-TC023 (Ø5.5, H=2mm)
CW-TC033 (Ø5.5, H=3mm)



CW-TCI13 (Ø5.5, H=1mm)
CW-TCI23 (Ø5.5, H=2mm)
CW-TCI33 (Ø5.5, H=3mm)



CW-TPH50 (Ø4.5mm)



CW-TPA10 (Ø5.5, H=1mm)

Gold Plastic Cylinders



CW-GPC10 (Ø5.5mm H=1mm)



CW-GP010 (Ø5.5mm H=1mm)

Abutments for Cement Retained Restorations



CW-SMAC1 (Ø5.5mm H=1mm)
CW-MC103 (Ø5.5mm H=3mm)



CW-AT001 (Ø5.5mm H=1mm)
CW-AT002 (Ø5.5mm H=2mm)
CW-AT003 (Ø5.5mm H=3mm)



CW-CPK61 C=6mm H=1mm
CW-CPK62 C=6mm H=2mm
CW-CPK63 C=6mm H=3mm
CW-CPK64 C=6mm H=4mm

CW-CPK81 C=8mm H=1mm
CW-CPK82 C=8mm H=2mm
CW-CPK83 C=8mm H=3mm
CW-CPK84 C=8mm H=4mm



MW-IC040 MW-ICH40

* Includes burnout anti-rotation plastic cap MW-ICH40

Angulated Abutment



CW-AN151 (Ø5.5mm H=1mm)
CW-AT152 (Ø5.5mm H=2mm)
CW-AT153 (Ø5.5mm H=3mm)



CW-AN251 (Ø5.5mm H=1mm)
CW-AT252 (Ø5.5mm H=2mm)
CW-AT253 (Ø5.5mm H=3mm)

Zirconia Solution



Non-Cemented

CK-WPZ4 H=0.5mm

* Available in kits only



Pre-Cemented

CK-WPC4 H=0.5mm

* Available in kits only

○ Free rotation

⊗ Anti-rotation

*All superstructures are provided with the relevant prosthetic screws

WIDE

Multi-Unit



CM-W1480 (H=1mm)
CM-W3480 (H=3mm)
CM-W5480 (H=5mm)

Angulated



CM-WA171
CM-WA301

Ball Attachments



CW-B1500 (H=1mm)
CW-B3500 (H=3mm)
CW-B5500 (H=5mm)

OT-Equators



CK-WOE1 (H=1mm)
CK-WOE3 (H=3mm)
CK-WOE5 (H=5mm)

* Available in kits only

Ti-Base



CW-TB001 (Ø5.7mm H=0.5mm)
CW-TB015 (Ø5.7mm H=1.5mm)



CW-TB002 (Ø5.7mm H=0.5mm)
CW-TBF15 (Ø5.7mm H=1.5mm)

Analog



CW-RSM10

Burnout



GI-WCH11


* Scanning Abutment



CW-SP001

* Scanning Abutment is sold separately and is not part of the Ti-Base package.

Implant Details



5.0
Diameter

*All superstructures are provided
with the relevant prosthetic screws

- Free rotation
- ◻ Anti-rotation




Healing Cap



MU-H4480

Impression Copings



MU-I0480

Analog



MU-RSM48

Temporary Cylinder



MU-T0480



MU-S0220

Plastic Cylinders



MU-CP048



MU-S0220

Gold Plastic Cylinders



MU-MG480



MU-G0220



MU-BGC48



Ball Attachments



MB-RS375



OT-Equators



OE-IC002



OE-RSM60

Prosthetic Options (Complementary Components)

Gold Ball Cap	Plastic Ball Cap	Titanium Ball Cap	Metal Housing & Plastic Caps	Metal Housing & Transfer Caps	Impression Transfer Cap	CAD/CAM Cementing Cap	Scanning Abutment
						 MU-CC480	 MU-SP001
 MB-GPS10	 MB-SF200	 MB-TTP10	 MB-DMH10  MB-PPC10  MB-PPR10				
				 OE-RCW01 OE-RCY01 OE-MH001 OE-RCB01 OE-RCP01 OE-RCV01	 OE-IC001		

Restorative Components

Guidelines for usage of MIS Restorative Components

Location ►	Anterior Maxilla	Canine, Premolars and Molars				Premolars and Molars
Crown/ Implant Inclination Ratio	15° or 25°	Up to 3.5°			15° or 25°	Up to 4°
Gingival Height	Up to 3mm buccal Up to 4mm palatal	Up to 3mm buccal Up to 4mm lingual/palatal	Adjustable, based on abutment preparation	Up to 2mm	Up to 3mm buccal Up to 4mm lingual	According to gingival height available in heights of 1,2,3,4mm
Catalog Number	CS-A1510 CS-A2510	CS-A0010	CS-MAC10	CS-GPC10 CS-GP010	CS-AN151 CS-AN251	CS-CPK61 CS-CPK62 CS-CPK63 CS-CPK64
Abutment ► description	Esthetic angulated abutment	Esthetic abutment	Standard post abutment	Screw-retained gold plastic abutment	Angulated abutment	Anatomic transgingival abutment
						
						

*Recommendation only

